es [INSTRUCCIONES DE USO]

Mini-Implantes dentales

Datos del fabricante:

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3. Intended use / Purpose:

Mini-implants are indicated to be fixed in the jawbone, periodontal tissue or palate. They are used as a temporary anchorage point in orthodontic/disjunction treatments.

4. Medical purpose / Indication:

Dental mini-implants are indicated as an anchor point in orthodontic treatments, intended to be used as a fixed anchor point for fixing devices that facilitate orthodontic movement of the teeth, such as the intrusion of molars, the retraction of incisors and the intrusion of canines, among others. Also, mini-implants are used in procedures that aim to expand the maxilla, separating the midpalatal suture.

Correct treatment that requires the use of orthodontic/disjunction mini-implants must be based on a diagnosis, planning, and treatment design that meet the needs and expectations of the patient.

5. Intended user:

The intended user is a medical professional, responsible for appropriate orthodontic surgical procedures and techniques; who must evaluate the suitability of the procedure used, based on their personal medical knowledge, education, training and experience (at least 1 year of experience in dental treatments) from the process of fabrication hasta el momento de su distribución. Este producto está destinado a ser utilizado en un ambiente estéril.

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The intended patient requires this type of treatment to prevent and/or correct the different alterations in dental development, the position of the jaws, the shapes of the dental arches and/or separate the cleft palate by expanding the jaw to increase the length of the top archade and to gain space; all with the aim of restoring the functional and morphological balance of the mouth and face. See the "Contraindications" section.

7. Indications for use:

It is recommended to follow the instructions for use recommended by the manufacturer, as well as to take the necessary protective and clothing measures. Improper use can lead to injure the patient.

To select the type of mini-implant according to its diameter, length, type of head and neck depending on its intended purpose or location.

Using the manual handle, attach the mini-implant into the insert.

To check the type of bone tissue, and the soft tissue to be implanted, and anaesthetize the area.

If necessary, to mark the position of the mini-implant in the soft tissue and make an incision.

Place the screw directly on the mucosa, and then press and rotate the manual handle, so that the tip of the screw can dig a hole in the cortical bone surface. (While inserting the screw, turn it slowly to prevent the insert and manual handle from vibrating).

Stop inserting the mini-implant when the lower end of the head of the mini-implant touches the mucosa and gently remove the insert and the manual handle.

The design of the mini-implants enables them to be used together with other products and accessories commonly used in orthodontics/disjunction treatments, such as wires, metal ligatures, elastic modules, chains, springs, etc.

This product is intended to be used in a sterile environment. Any serious incident related to the device must be notified to the manufacturer and to the competent authority of the Member State in which the intended user and/or the patient is established.

Before use, check that there is no deformation or damage to the product.

If during implantation is detected a variation in the consistency of the bone, stop the insertion immediately and make sure that the mini-implant is not inserted into unsuitable bone or dental structures such as the root of the tooth.

Excessive perpendicularly force on the alveolar bone can cause the mini-implant to tilt and can lead to its fracture. The maximum insertion depth is up to the neck of the mini-implant. Pay special attention not to apply excessive force or torque during insertion.

To prevent infections around the implantable area, it must be kept clean.

8. Contraindications:

It is contraindicated to use the mini-implants in patients: Allergic or hypersensitive to the materials used.

With severe or uncontrolled periodontal diseases.

With limited bone availability and density, which prevents the retention of the mini-implant along the treatment.

Who suffer from osteoporosis, have inhibited revascularization and/or deficient bone formation.

Inadequate oral habits, with infection, diabetics and/or smokers.

Because they can cause this medical device to loosen, bend and break, or not be fixed properly in the bone.

9. Possible adverse reactions and risks:

An improper use of the products or an inadequate implantation protocol can cause surgical or postoperative pain, inflammation/infection of the mucosa in the area, bruising or bleeding, mobility of the device in the bone, loss or fall of the mini-implant, alveolar fracture, mini-implant fracture, damage to adjacent roots or teeth, injury to adjacent anatomical structures (nasal cavities, maxillary sinus, vascular/nerve structures).

Recommended ideal storage conditions: mini-implants should be stored in its original package in a clean and dry storage areas, not exposed to direct sunlight for 5 years (expiration period to ensure its sterility).

Incorrect storage during transport can influence the characteristics of the device and lead to its failure.

2. Description and material:

Mini-implants are medical devices (dental), made of commercially pure Titanium Grade IV cold drawn bars. Orthodontic/disjunction mini-implants are medical devices class IIb.

For more information, please visit our website. [*]

3. Intended use / Purpose:

Mini-implants are indicated to be fixed in the jawbone, periodontal tissue or palate. They are used as a temporary anchorage point in orthodontic/disjunction treatments.

4. Medical purpose / Indication:

Dental mini-implants are indicated as an anchor point in orthodontic treatments, intended to be used as a fixed anchor point for fixing devices that facilitate orthodontic movement of the teeth, such as the intrusion of molars, the retraction of incisors and the intrusion of canines, among others.

Also, mini-implants are used in procedures that aim to expand the maxilla, separating the midpalatal suture.

Correct treatment that requires the use of orthodontic/disjunction mini-implants must be based on a diagnosis, planning, and treatment design that meet the needs and expectations of the patient.

5. Intended user:

The intended user is a medical professional, responsible for appropriate orthodontic surgical procedures and techniques; who must evaluate the suitability of the procedure used, based on their personal medical knowledge, education, training and experience (at least 1 year of experience in dental treatments) from the process of fabrication hasta el momento de su distribución.

This product is indicated to be sterilized (the irradiation label must be red).

It is recommended that the intended user has received regulated and supervised training by expert professionals on orthodontic/disjunction treatment techniques; the intended user should receive personalized assistance and the intended user should attend special training before using the products. Product specifications alone do not guarantee proper use. MPI offers courses for different levels of knowledge and experience.

6. Intended patient:

The intended user is a medical professional, responsible for appropriate orthodontic surgical procedures and techniques; who must evaluate the suitability of the procedure used, based on their personal medical knowledge, education, training and experience for use.

The intended user who uses these products must have experience to use them safely and properly according to these instructions for use, in addition the intended user must inform intended patients about contraindications, risk and possible adverse reactions, as well as the importance of maintaining good dental hygiene.

The intended patient requires this type of treatment to prevent and/or correct the different alterations in dental development, the case, the indications and contraindications of these instructions for use will be considered in order to minimize error variables.

Once the complete case has been planned, it is necessary to study the anatomical structure of the area where the mini-implants are going to be placed, and decide the type, diameter and length, to minimize possible surgical damage.

Always, take into account the sterilization conditions of the operative field.

7. Description and material:

Mini-implants are medical devices (dental) manufactured in titanium with commercial purity of Grade IV, cold drawn. Mini-implants of orthodontic/disjunction are medical devices class IIb.

Its sterilization is carried out in accordance with the instructions for use.

It is recommended that the intended user has received training and supervision by expert professionals on orthodontic/disjunction treatment techniques; the intended user should receive personalized assistance and the intended user should attend special training before using the products. Product specifications alone do not guarantee proper use. MPI offers courses for different levels of knowledge and experience.

The intended user who uses these products must have experience to use them safely and properly according to these instructions for use, in addition the intended user must inform intended patients about contraindications, risk and possible adverse reactions, as well as the importance of maintaining good dental hygiene.

The intended patient requires this type of treatment to prevent and/or correct the different alterations in dental development, the case, the indications and contraindications of these instructions for use will be considered in order to minimize error variables.

Once the complete case has been planned, it is necessary to study the anatomical structure of the area where the mini-implants are going to be placed, and decide the type, diameter and length, to minimize possible surgical damage.

Always, take into account the sterilization conditions of the operative field.

When performing surgery, always consider the anatomical structures adjacent to the mini-implant, such as the nasal cavity, maxillary sinus, vascular/nerve structures.

The intended user who uses these products must have experience to use them safely and properly according to these instructions for use, in addition the intended user must inform intended patients about contraindications, risk and possible adverse reactions, as well as the importance of maintaining good dental hygiene.

The intended patient requires this type of treatment to prevent and/or correct the different alterations in dental development, the case, the indications and contraindications of these instructions for use will be considered in order to minimize error variables.

Once the complete case has been planned, it is necessary to study the anatomical structure of the area where the mini-implants are going to be placed, and decide the type, diameter and length, to minimize possible surgical damage.

Always, take into account the sterilization conditions of the operative field.

When performing surgery, always consider the anatomical structures adjacent to the mini-implant, such as the nasal cavity, maxillary sinus, vascular/nerve structures.

The intended user who uses these products must have experience to use them safely and properly according to these instructions for use, in addition the intended user must inform intended patients about contraindications, risk and possible adverse reactions, as well as the importance of maintaining good dental hygiene.

The intended patient requires this type of treatment to prevent and/or correct the different alterations in dental development, the case, the indications and contraindications of these instructions for use will be considered in order to minimize error variables.

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