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SIEWERT BRIDGE™ SOLUTIONS
C/ Aquilón, 2
28223 Pozuelo de Alarcón
Madrid, España

(+34) 91 711 42 65
www.siewertbridge.com

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
1. Product:

Dental Prosthetic components (Abutment Screw, Laboratory screw, Scan Body), and dental instrument (Screwdriver laboratory) to use with the SIEWERT BRIDGE™ Solutions.

2. Description and material:

They are dental devices designed to be used as a support in prosthetic restoration.

Product type	Abutment screw (SIEWERT BRIDGE™ SCREW Multi Unit Patient)
Reference	SBSMUP
Description	Designed for fixing the patient's prosthesis. The maximum torque recommended is 15 N-cm. Use: Definitive. It is shipped sterile packaged in a blister. It is a medical device class IIb.
Material	Titanium Grade V with DLC coating.

Product type	Abutment screw (SIEWERT BRIDGE™ SCREW Multi Unit Patient Angled)
Reference	SBSMUPA
Description	Designed for fixing the patient's prosthesis. The maximum torque recommended is 15 N-cm. Use: Definitive. It is shipped sterile packaged in a blister. It is a medical device class IIb.
Material	Titanium Grade V with DLC coating.

Product type	Laboratory screw (SIEWERT BRIDGE™ SCREW Multi Unit Laboratory)
Reference	SBSMUL
Description	Designed for the lab process and for the fixation of the scan bodies. It is not a medical device. It is shipped packaged in a plastic pouch, non-sterile, and it can be reused. It is used only in a dental laboratory (no intraoral use) and it has not cleaning and sterilization requirements.
Material	Titanium Grade V.

Product type	Scan Body (SIEWERT BRIDGE™ Multi Unit Scan Body)
Reference	SBMUUSB
Description	It is a device designed to establish the virtual 3D position of the implant in the oral cavity. Use: Temporary. It is a medical device class I. It is shipped packaged in a plastic pouch, non-sterile, and it can be reused. It is used only in a dental laboratory (no intraoral use) and it has not cleaning and sterilization requirements.
Material	PEEK.

Product type	ScrewDriver Laboratory
Referencia	SDLB
Description	It is a reusable dental instrument. It is used for fixing the pieces. It is not a medical device. It is shipped packaged in a plastic pouch, non-sterile, and it can be reused. It is used only in a dental laboratory (no intraoral use) and it has not cleaning and sterilization requirements.
Material	Stainless Steel.

3. Intended use / Purpose:

The SIEWERT BRIDGE™ Solutions allows extensive control by creating an exact PMMA model of the final denture first for temporary trials. This way possible errors are identified before the final and definite denture is manufactured. The SIEWERT BRIDGE™ can therefore be used on practically all implant systems.

The definitive SIEWERT BRIDGE™ avoids screw loosening or screw fracture; it is used to avoid chipping in the ceramic interface; to achieve stable cervical bone around the implants; it avoids changes due to abrasion in the occlusion; and it enables a passive fit upon the multi-unit abutments in the patient's mouth.

4. Medical purpose / Indication:

The SIEWERT BRIDGE™ Solutions is indicated for partial or total rehabilitation of edentulous spans using an immediate or delayed loading protocol, according to the intended user's choice, the technique applied and the primary stability achieved.

5. Intended user:

The intended user is a medical professional, responsible for the appropriate surgical procedures and restorative techniques; who must evaluate the suitability of the procedure used, based on their personal medical knowledge, education, training and experience (at least 1 year of experience in dental surgeries with this type of product).

6. Intended patient:

Adult patients with unitary spaces without tooth with adjacent and healthy teeth and/or partial edentulism with missing posterior tooth and/or complete edentulism; that cannot tolerate removable prostheses, with high aesthetic and/or functional demands, and who are not contraindicated. See chapter "Contraindications".

7. Indications for use:

-It is recommended to follow the instructions for use recommended by the manufacturer, as well as to take the necessary protective and clothing measures. Improper use can lead to injure the patient.

-Intended users should review the expiration date indicated on the package and check that the package is not damaged or opened to maintain the sterility.

-Before placing the piece, the intended user must verify that it is compatible with the implant model placed and the appropriate conditions for using them.

-When it is decided to discard a piece, the intended user must dispose it safely following the protocol for the elimination of medical devices defined in the dental clinic, in accordance with hospital regulations and local existing legislation, to guarantee the safety elimination and avoid infections and contamination. [*]

8. Contraindications:

Contraindicated patients:

-Not suitable for an oral surgical procedure (for example, patients who have undergone radiotherapy or who have been treated with bisphosphonates, should contact the specialist who has treated them to verify the risks and contraindications). -Allergic or hypersensitive to the materials used.

-With periodontal diseases, limited bone availability and density, diabetics and/or smokers and/or bruxism. In case of a patient with bruxism (the involuntary habit of clenching and grinding the teeth without functional purpose), the intended user should assess it in the diagnosis and also in the treatment plan because it should be reflected in the occlusal protection of the implant during bone integration, and in the design of the final prosthesis. -Patients who cannot undergo surgery for any medical condition. None have been indicated for laboratory components.

The intended user must inform the intended patient of these contraindications.

9. Possible adverse reactions and risks:

-The inadequate application of higher torque values than the recommended ones, could damage the connection between the prosthetic component and the implant; in addition, it could lead to a lack of primary stability.

-The unskilled or untrained user, without planning clinical cases and/or in contraindicated patients, may cause harm to the patient.

-An inappropriate environment where surgery is to be performed/ improper storage conditions, can cause contamination of the product,

and can lead to infection of the patient and/or non-osseointegration of the implant.

To date, no other possible significant adverse effects caused by these products are known.

The intended user must inform the intended patient about the importance of maintaining good dental hygiene, and about these possible adverse reactions and risks.

10. Sterilization:

Some of the indicated products are delivered sterile and others are delivered non-sterile (for use in dental laboratory only, no intraoral use).

Abutment screw	Supplied sterile: Sterilized by gamma irradiation
Laboratory screw	Supplied Non-sterile
Scan Body	Supplied Non-sterile
ScrewDriver Laboratory	Supplied Non-sterile

11. Precautions:

-Before placing the product, check if it must be sterile according to the "Sterilization" section; in case it must be sterile check that it has been sterilized (the irradiation label must be red).

-It is recommended that the intended user have received regulated and supervised training by expert professionals on surgical techniques in oral implantology; the intended user should consult the specific surgical and prosthetic protocols; the intended user should attend special training before using the products. MPI offers courses for different levels of knowledge and experience.

-The intended user who uses these products must have experience to use them safely and properly according to these instructions for use, in addition the intended user must inform intended patients about contraindications and possible adverse reactions, as well as the importance of maintaining good dental hygiene.

-During the phase of the diagnosis and planning the surgery, the indications and contraindications of this instructions for use will be considered, the MPI implant system and its prosthetic possibilities should be known, in order to minimize error variables. [*]

-It is recommended to follow the instructions for use recommended by the manufacturer, as well as to take the necessary protective and clothing measures. Improper use can lead to injure the patient.

-Intended users should review the expiration date indicated on the package and check that the package is not damaged or opened to maintain the sterility.

-Before placing the piece, the intended user must verify that it is compatible with the implant model placed and the appropriate conditions for using them.

-When it is decided to discard a piece, the intended user must dispose it safely following the protocol for the elimination of medical devices defined in the dental clinic, in accordance with hospital regulations and local existing legislation, to guarantee the safety elimination and avoid infections and contamination. [*]

-If the products (indicated as single use devices) are re-used/ re-sterilized, the sterilization levels offered by MPI are not guaranteed, as they can cause infections and adverse reactions in the patient and therefore the non-osseointegration of the implant.

-The product is perfectly packaged and heat-sealed. Check if the sealed package is damaged or opened because the sterility conditions of the product (in case of sterile products) and cleaning conditions of the product could have been altered; in that case, do not use it.

-All products have a batch number to allow their traceability, providing greater security to the intended user, due to with this batch number it is possible to know the entire history of the product from the manufacturing process to the moment of its distribution.

-Any serious incident related to the device must be notified to the manufacturer and to the competent authority of the Member State in which the intended user and/or the patient is established.

13. Handling:

To place the prosthetic components, use compatible dental instruments and follow the recommendations in these instructions to ensure compatibility with the SIEWERT BRIDGE™ Solutions.

14. Storage and transport:

The products are wrapped and packaged according to defined and validated protocols that guarantee optimal results. The original packaging has been validated simulating the transport and storage conditions to guarantee the cleanliness of the packaged products, their protection against impacts, and the sterility of them (in case of sterile products).

Incorrect storage during transport can influence the characteristics of the device and lead to its failure.

Recommended ideal storage conditions:

-Sterile products should be stored in its original package in a clean and dry storage areas, not exposed to direct sunlight for 5 years (expiration period to ensure its sterility).

-Non-sterile products should be stored in its original package in a clean and dry storage areas, not exposed to direct sunlight.

15. Labelling symbols:

	Manufacturer.
	Use-by date.
	Batch code.
	Catalogue number.
	País de fabricación con el código de dos letras del país definido en la Norma ISO 3166-1. La fecha de fabricación se puede añadir adyacente al símbolo.
	Sterilized using irradiation.
	Do not re-sterilize.
	Do not use if package is damaged and consult instructions for use.
	Keep away from sunlight.
	Keep dry.
	Unique device identification.
	Do not re-use.
	Consult instructions for use.
	Caution.
	Medical device.
	CE symbol with the notified body number 1639.
	CE symbol.
	Unique Device Identification.
	Non-sterile.
	UDI - DataMatrix with the information: (01) UDI-DI (17) Expiration date (YYMMDD) (10) Número de lote

[*] For more information, please consult the catalogue and the following websites:
<https://www.mpimplants.com>, <https://www.siewertbridge.com>.

Additional Information: See the "Summary of Safety and Clinical Performance (SSCP)" at www.mpimplants.com.

Updated revision: April 2023 – Last revision: 4