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en INSTRUCTIONS FOR USE / CLASS IIB PROSTHETIC COMPONENTS

1. Product:
Dental prosthetic components MPI Privilege® External Connection, Dental prosthetic components MPI Privilege® Internal Connection, Dental prosthetic components MPI Excellence® Internal Connection, Dental prosthetic components MPI Short® Internal Connection and Dental prosthetic components MPI All-In® Internal Connection.

Models:
Abutment screw, Cover screw, Healing abutment, Transepithelial abutment, Conical abutment, Transepithelial screw, Overdenture abutment, Ti-Base, Cr-Co Base castable abutment.

2. Description and material:

They are medical devices (dental), intended for single use only, with direct connection to the abutment or dental structure, designed to be used as a support in prosthetic restoration.
They are medical devices class IIB. See table attached at the end of the instructions.
The Grade V Titanium of the dental prosthetic components manufactured by MPI meets the chemical composition requirements detailed in ISO 5832-3. See table attached at the end of the instructions.
The CoCr of the dental prosthetic components manufactured by MPI meets the chemical composition requirements detailed in ISO 5832-12. See table attached at the end of the instructions.

3. Intended purpose:

Our dental prosthetic components are intended for single use only, their intended purpose is to create a structure on a dental implant in the maxilla or jaw, in order to restore masticatory function, to achieve an aesthetic and functional improvement, and a decrease in bone loss.

4. Medical purpose / Indication:

Dental prosthetic components are indicated for partial or total rehabilitation of edentulous spans using an immediate or delayed loading protocol, according to the intended user's choice, the technique applied and the primary stability achieved.

5. Intended user:

The intended user is a medical professional, responsible for the appropriate surgical procedures and restorative techniques; who must evaluate the suitability of the procedure used, based on their personal medical knowledge, education, training and experience (at least 1 year of experience in dental surgeries with this type of product).

6. Intended patient:

Dental implants and prosthetic components are indicated to be placed in:
Adult patients (generally older than 17-18 years) with (jaw)maxilla growth complete, with unitary spaces without tooth with adjacent and healthy teeth and/or partial edentulism with missing posterior tooth and/or complete edentulism; that cannot tolerate removable prosthesis, with high aesthetic and/or functional demands, and who are not contraindicated.

See chapter "Contraindications".

There is no available data related to special patient populations such as children and during pregnancy or breastfeeding. Therefore, the use of MPI Prosthetic Components must be according to the criteria of the specialist healthcare professional.

7. Indications for use:

- It is recommended to follow the instructions for use recommended by the manufacturer, as well as to take the necessary protective and clothing measures. Improper use can lead to injure the patient.
- They are intended to be used in a sterile environment.
- The prosthetic components indicated must be used with the rest of the articles and references of the system to which they belong in order to guarantee a good result.
- Intended users should review the expiration date indicated on the package and check that the package is not damaged or opened to maintain the cleaning of the packaged product, and the sterility (in case of applicable).
- These products are available in different lengths and diameters in order to be able to rehabilitate different cases.
- Before placing the prosthetic component, the intended user must verify that it is compatible with the implant model placed and the appropriate conditions for using them.
- When it is decided to discard a piece, the intended user must dispose it safely following the protocol for the elimination of medical devices defined in the dental clinic, in accordance with hospital rules and local existing legislation, to guarantee the safety elimination and avoid infections and contamination.
- The maximum torques recommended for each type of component (established according to the mechanical tests carried out), are indicated in the "Description and material" section.

8. Contraindications:

- It is contraindicated to use the dental implants and prosthetic components in patients:
 - Not suitable for an oral surgical procedure (for example, patients who have undergone radiotherapy or who have been treated with bisphosphonates, should contact the specialist who has treated them to verify the risks and contraindications).
 - Allergic or hypersensitive to the materials used.
 - With periodontal diseases, limited bone availability and density, diabetics and/or smokers and/or bruxism. In case of a patient with bruxism (the involuntary habit of clenching and grinding the teeth without functional purpose), the intended user should assess it in the diagnosis and also in the treatment plan because it should be reflected in the occlusal protection of the implant during bone integration, and in the design of the final prosthesis.
- The placement of dental implants and prosthetic components is contraindicated in those patients who cannot undergo surgery for any medical condition.
- The intended user must inform the intended patient of these contraindications.

9. Possible adverse reactions and risks:

- An improper use of the products or an inadequate surgical technique, can cause bone loss, peri-implant disease, inflammation, discomfort and pain in the patient, and the consequent non-oseointegration of the implant.
- The inadequate application of torque values could damage the connection between the dental implant and prosthetic component causing abutment breakage or screw loosening; in addition, it can lead to a lack of primary stability, which would mean failure in the implant osseointegration.
- The unskilled or untrained user, without planning clinical cases and/or in contraindicated patients, may cause harm to the patient.
- An inappropriate environment where surgery is to be performed/improper storage conditions, can cause contamination of the product, and can lead to infection of the patient and/or non-oseointegration of the implant.
- To date, no other possible significant adverse effects caused by these products are known.
- The intended user must inform the intended patient about the importance of maintaining good dental hygiene, and about these possible adverse reactions and risks.

10. Sterilization:

Some of the indicated products are delivered sterile (those that will be used in a dental clinic), and others are delivered non-sterile (those that will be used in a dental laboratory). See table with description at the end of the instructions.

11. Precautions:

- Before placing the product, check if it must be sterile according to the "Sterilization" section; in case it must be sterile check that it has been sterilized (the irradiation label must be red).
- It is recommended that the intended user has received regulated and supervised training by expert professionals on surgical techniques in oral implantology; the intended user should consult the specific surgical and prosthetic protocols; the intended user should receive personalized assistance and the intended user should attend special training before using the product. MPI offers courses for different levels of knowledge and experience.
- The intended user who uses these products must have experience to use them safely and properly according to these instructions for use, in addition the intended user must inform intended patients about contraindications and possible adverse reactions, as well as the importance of maintaining good dental hygiene.
- During the phase of diagnosis and planning the surgery, the indications and contraindications of this instructions for use will be considered, the MPI implant system and its prosthetic possibilities should be known, in order to minimize error variables.
- Once the complete case has been planned, and the location, the diameter and the length of the implant have been decided, the type of prosthesis to be used will be designed.
- Subsequently, the patient will undergo surgery, taking into account the design of the incision, the drilling protocol, the placement of the implant and the sterilization conditions of the operative field.

12. Warnings:

- It is necessary to consider the medical situation of the patient, the physiological and anatomical conditions, and the bone quality and bone quantity in the area where the surgery will be performed. An inadequate technique can incur any of the possible adverse reactions and risks mentioned above.
- This product is only for single use, it cannot be re-sterilized or re-used, if the products are re-used, the sterilization levels offered by MPI are not guaranteed, as they can cause infections and adverse reactions in the patient and therefore the non-oseointegration of the implant.
- The product is perfectly packaged and heat-sealed. Check if the sealed package is damaged or opened because the sterility conditions of the product (in case of sterile products) and cleaning conditions of the product could have been altered; in that case, do not use it.
- The product Cr-Co base castable abutment contains cobalt, a substance classified as CMR (carcinogenic, mutagenic, or toxic to reproduction). Its use should be carefully evaluated in sensitive patients or vulnerable groups.
- This product is intended to be used in a sterile environment.
- Any serious incident related to the device must be notified to the manufacturer and to the competent authority of the Member State in which the intended user and/or the patient is established.

13. Magnetic resonance safety information:

MPI Dental Prosthetic Components have not been evaluated for safety and compatibility in the magnetic resonance environment. Scanning a patient with this medical device may not be safe.

14. Handling:

- It is important to verify the compatibility of the implant model to be used.
- To place the prosthetic components, compatible dental instruments provided by MPI must be used; to activate the rotation of the power instruments, they must be connected to the handpiece (contra-angle) of the implant motor.

15. Storage and transport:

- The prosthetic components are wrapped and packaged according to defined and validated protocols that guarantee optimal results. The original packaging has been validated simulating the transport and storage conditions to guarantee the cleanliness of the packaged products, their protection against impacts, and the sterility of them (in case of sterile products).
- Recommended ideal storage conditions:
 - For sterile products: store in its original packaging in a clean and dry storage areas, not exposed to direct sunlight for 5 years (expiration period to ensure its sterility).
 - Non-sterile products should be stored in its original packaging in a clean and dry storage areas, not exposed to direct sunlight.
- Incorrect storage during transport can influence the characteristics of the device and lead to its failure.

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