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en

INSTRUCTIONS FOR USE / I CLASS SUPPORT FOR DENTAL PRODUCTS

Manufacturer data:
Name: MEDICAL PRECISION IMPLANTS S.A. , Address: c/ Empleo, 21, Pl. Los Olivos, 28906, Getafe, Madrid, Spain
Telephone: +34 91 684 60 63
Website: www.mpimplants.com

1. Product:
Support for dental products. **Types:** See table at the end of the instructions.

2. Description and material:
Container for storing dental products.
It is a medical device Class I in accordance with the Rule 1 - Annex VIII of the REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL.

3. Intended use/purpose:
It is designed to hold unwrapped dental components and instruments, during sterilization and for their subsequent storage.
It is intended to be used in dental implant placement procedures.
It is permeable to allow moisture to escape from within after the sterilization process is finished, whilst protecting the contents from being contaminated.
This is a reusable device.

4. Medical purpose/indication:
Placement of a dental implant in edentulous patient.

5. Intended user:
The intended user is a medical professional, responsible for the appropriate surgical procedures and restorative techniques, who should assess the suitability of the procedure used, based on their personal medical knowledge, education, training and experience (at least 1 year of documented experience in dental surgeries with these devices).

6. Intended patient:
Adult patients (in general 17-18 years and older) with jaw/maxilla growth complete, to be treated with a dental implant system and who are not contraindicated. See the chapter "contraindications".

7. Indications for use:
- It is recommended to follow the instructions for use recommended by the manufacturer, as well as to adopt the necessary protective and clothing measures. Inadequate use of the instruments may result in patient damage.
- Before using the instruments, the intended user must inspect the device to ensure that the integrity of the device is maintained, once verified their good condition the intended user must disinfect, clean and sterilize them by autoclave moist heat. The intended user is responsible for the inspection, disinfection, cleaning and sterilization of them.
- When deciding to replace the instruments, the intended users should dispose them following the protocol for disposal of medical devices defined in the dental clinic, to ensure their safe disposal.

8. Contraindications:
It is contraindicated to use the products in patients:
- Not suitable for an oral surgical procedure (for example, patients who have undergone radiotherapy or who have been treated with bisphosphonates, should contact the specialist who has treated them to verify the risks and contraindications).
- Allergic or hypersensitive to the materials used.
- Contraindicated for treatment with dental implants or restorative components.
- With periodontal diseases, limited bone availability and density, diabetes and / or smokers.
The intended user must inform the intended patient about these contraindications.

9. Possible adverse reactions and risks:
- Mishandling/inadequate use of devices can cause the injury of dimensional and ergonomic features of the instrument.
- Unqualified or untrained user who use the devices without planning clinical cases and/or in contraindicated patients can cause patient harm.
- Inappropriate environment in which the devices are intended to be used, an inappropriate/incorrect disinfection, cleaning and sterilization of devices and incorrect storage conditions can cause device contamination, infection of the patient, and/or wear effect that can cause nerve damage.
To date, no other possible significant adverse effects caused by these products are known.
The intended user must inform the intended patient about the importance of maintaining good dental hygiene, and about these possible adverse reactions/risks.

10. Cleaning, disinfection and sterilization:
- The devices are delivered non-sterile and they are intended to be reused.
- The cleaning, disinfection and sterilization stages must be respected before the first use and each time the devices are reused. The intended user is responsible for its correct cleaning, disinfection and sterilization.
- Before sterilization by autoclave moist heat, to guarantee the elimination of contaminating microorganisms, the devices must be disinfected and cleaned with cleaning and disinfection agents with CE certificate, with bactericidal, fungicidal and virucidal properties, according to the instructions manufacturer.

Step 1: Pre-disinfection
- After using the devices in the surgery, to eliminate the residues and remains of blood, and to prevent that the blood will be dried on its surface and thus facilitate their subsequent cleaning, immerse them in a disinfectant solution (with EC certificate, with bactericidal, fungicidal and virucidal properties, with any substance not compatible with the drills), according to the manufacturer's instructions (pay attention to concentrations and immersion time; excessive concentration may cause corrosion of the products).
- Disinfectant solutions containing substances not compatible with the devices should not be used.
- Rinse the devices with purified water. If they have visible impurities, it is recommended to brush them manually.

Step 2: Manual cleaning and disinfection
- Immerse the devices in a cleaning and disinfectant solution (with EC certificate, with bactericidal, fungicidal and virucidal properties), in an ultrasonic bath, according to the manufacturer's instructions.
- Immerse the devices in an ultrasonic bath for 15 min at 35°C, maintaining a maximum temperature of 45°C to avoid protein coagulation.
- Repeat the process again with new purified water, to make a total of two rinses.
- Remove the devices and rinse them (immerse them, shake them and keep them for 1 minute in purified water).
- Dry the products with filtered air.

Step 3: Inspection
- Inspect devices and discard those that have defects that may affect their resistance, safety or performance.
- Wash and disinfect dirty devices again.

Step 4: Sterilization
- To ensure the complete elimination of all life forms, including spores, the intended user must perform the sterilization of the devices before using them.
- To sterilize the devices, store them in a suitable sterilization package, seal it, and sterilize them in autoclave moist heat at 134°C for 18 minutes.
- To ensure that the sterilizer maintains adequate working conditions, it is important to inspect and clean it regularly. It is recommended periodic tests, cleaning and calibration following the instructions for use and maintenance of the sterilizer.

11. Precautions:
- The intended user who uses the product must have experience to use it safely and adequately to these instructions for use.
- It is recommended that the intended user consult the surgical protocols, receive personalized assistance, and attend special training before using the devices. MPI offers courses for different levels of knowledge and experience.
For more information, please see our website. [*]

12. Warnings:
- After each use, clean and disinfect the devices. Improper cleaning can lead to improper sterilization. The use of hydrogen peroxide or other oxidizing agents will damage the surface of the devices. Use disinfectant solutions approved for their effectiveness.
- Dry heat sterilization is not allowed, as it directly affects the performance of the devices.
- Ensure that the devices are completely dry. If they do not dry completely during sterilization, they can remain wet and cause discoloration and oxidation. The oxidized instruments can contaminate the water circuit, and the oxide particles can cause the initial oxidation of the initial instruments in all future sterilization cycles.
- Do not use the instruments if the sealed package in which the instrument is shipped is damaged or opened.
- These devices are intended to be used in a sterile environment.
- Any serious incident related to the device must be notified to the manufacturer and to the competent authority of the Member State in which the intended user and / or the patient is established.
- For best results, these instruments should only be used with the indicated implant systems.

13. Handling:
During the surgical procedure, the box should be placed on a sterile surface and handled with a sterile glove. Remove the cover to access the devices. The container itself indicates where to place each product.

The dental case "DC68" is designed to be used with the insert:
- "IDCP68", for the instruments to be used with the family MPI Privilege®.
- "IDCC68", for the instruments to be used with the families MPI Privilege® y MPI Excellence®.
In these, once the lid of the container is open, the tray can be removed. In the lower part there are three holes to place more instruments.

The dental case "DC36" is designed to be used with the insert:
- "IDC36", for the instruments to be used with all the MPI families.
- "IDCP36", for the instruments to be used with the family MPI Privilege®.
- "IDCE36", for the instruments to be used with the family MPI Excellence®.
- "IDCS36", for the instruments to be used with the family MPI Short®.
- "IDCA36", for the instruments to be used with the family MPI All-In®.

The dental case "STD P CASE" is designed to be used with the drill stoppers of the family MPI Privilege®.

The dental case "STDC CASE" is designed to be used with the drill stoppers of the conical drills of the families MPI Excellence®, MPI Short®, MPI All-In®.
To fit the drill with the lid stopper, press down on the corresponding hole in the case.

The dental case "DC001" is designed to be used with the insert "IDC001".
The container "CCT" is designed to be used with impactation pins and bone graft screws.

The container "POT-EXT" with the insert "BLOC-EX", is designed to be used with the extractors.

14. Storage and transport:
The original packaging has been validated simulating transport and storage conditions to guarantee the cleanliness of the devices packaged, and its protection against impacts.
After sterilization, keep the instruments in the sterilization packaging, storing them in an environment not exposed to direct sunlight.
Incorrect storage during transport can influence the characteristics of the device and lead to its failure.

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