**es** - **INSTRUCCIONES DE USO** / **en** - **INSTRUCTIONS FOR USE** / **fr** - **MODE D'EMPLOI**

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**en** **INSTRUCTIONS FOR USE** / **CLASS I PROSTHETIC COMPONENTS**

1.**Product:** Prosthetic components to use with the MPI Privilege® Implant system: external and internal connection, MPI Excellence®, MPI Short® and MPI

All-In®.

**Types:** Impression coping, Temporary abutment, Scanbody. [\*]

2.**Description and material:** They are medical devices (dental), with direct connection to the implant or dental structure, designed to be used as a support in prosthetic restoration. They are medical devices class I.

**Product type:** Impression coping **Description:** It is screwed directly to the implant when it has been integrated, for impression taking. **Use:** Transient, less than 60 minutes. **Types:** External and internal connection. **Open tray (C / A):** there is a short stem to be used as a C / C. **Closed tray (C / C):** the post has 3 flat sides to facilitate its positioning in the silicone, since it is a drag element. It is shipped packed in a pouch, not sterile. **Material:** Stainless steel.

**Product type:** Temporary abutment **Description:** It is screwed directly to the implant to make a provisional prosthesis. **Use:** Short term, less than 30 days. **Types:** External and internal connection. **Hexagonal:** Non-rotating and Cylindrical: Non-rotating. **Recommendations:** The biocompatible material (PEEK) is not recommended for casting. It is shipped packed in a pouch, not sterile. **Material:** PEEK.

**Product type:** Scanbody **Description:** It is a device designed to establish the virtual 3D position of the implant in relation to the oral cavity. **Use:** Transient, less than 60 minutes. It is shipped packed in a pouch, not sterile. **Material:** Titanium grade 5.

3. **Intended use / Purpose:** Dental prosthetic components indicated are designed and conceived to create a structure on a dental implant in the maxilla or jaw, in order to restore masticatory function, to achieve an aesthetic and functional improvement, and a decrease in bone loss.

4. **Medical purpose / Indication:** Dental prosthetic components are indicated for partial or total rehabilitation of edentulous spans using an immediate or delayed loading protocol, according to the intended user’s choice, the technique applied and the primary stability achieved.

5. **Intended user:** The intended user is a medical professional, responsible for the appropriate surgical procedures and restorative techniques; who must evaluate the suitability of the procedure used, based on their personal medical knowledge, education, training and experience (at least 1 year of experience in dental surgeries with this type of product).

6. **Intended patient:** Dental implants and prosthetic components are indicated to be placed in: Adult patients (generally older than 17-18 years) with jaw/maxilla growth complete, with unitary spaces without tooth with adjacent and healthy teeth and/or partial edentulism with missing posterior tooth and/or complete edentulism; that cannot tolerate removable prosthesis, with high aesthetic and/or functional demands, and who are not contraindicated. See chapter “Contraindications”.

7. **Indications for use:**

- It is recommended to follow the instructions for use recommended by the manufacturer, as well as to take the necessary protective and clothing measures. Improper use can lead to injure the patient.
- The prosthetic patients indicated must be used with the rest of the articles and references of the system to which they belong in order to guarantee a good result.
- Intended users should check that the package is not damaged or opened to maintain the cleaning of the packaged product.
- These products are available in different lengths and diameters in order to be able to rehabilitate different cases.
- Before placing the prosthetic component, the intended user must verify that it is compatible with the implant model placed and the appropriate conditions for using them.
- When it is decided to discard a piece, the intended user must dispose it safely following the protocol for the elimination of medical devices defined in the dental clinic/laboratory, in accordance with regulations and local existing legislation, to guarantee the safety elimination and avoid infections and contamination.[\*]

8. **Contraindications:** It is contraindicated to use the prosthetic components in patients:

- Not suitable for an oral surgical procedure (for example, patients who have undergone radiotherapy or who have been treated with bisphosphonates, should contact the specialist who has treated them to verify the risks and contraindications).
- Allergic or hypersensitive to the materials used.
- With periodontal diseases, limited bone availability and density, diabetics and/or smokers and/or bruxism. In case of a patient with bruxism (the involuntary habit of clenching and grinding the teeth without functional purpose), the intended user should assess it in the diagnosis and also in the treatment plan because it should be reflected in the occlusal protection of the implant during bone integration, and in the design of the final prosthesis.
- The placement of prosthetic components is contraindicated in those patients who cannot undergo surgery for any medical condition.
- The intended user must inform the intended patient of these contraindications.

These contraindications are not applicable for those prosthetic components intended to be used only in the dental laboratory.

9. **Possible adverse reactions and risks:**

- The unskilled or untrained user, without planning clinical cases and/or in contraindicated patients, may cause harm to the patient.
- An inappropriate environment where surgery is to be performed/improper storage conditions, can cause contamination of the product, and can lead to infection of the patient and/or non-osseointegration of the implant.

To date, no other possible significant adverse effects caused by these products are known.
The intended user must inform the intended patient about the importance of maintaining good dental hygiene, and about these possible adverse reactions and risks.

10. **Precautions:**

- It is recommended that the intended user has received regulated and supervised training by expert professionals on surgical techniques in oral implantology; the intended user should consult the specific surgical and prosthetic protocols; the intended user should receive personalized assistance and the intended user should attend special training before using the products. MPI offers courses for different levels of knowledge and experience.
- The intended user who uses these products must have experience to use them safely and properly according to these instructions for use, in addition the intended user must inform intended patients about contraindications and possible adverse reactions, as well as the importance of maintaining good dental hygiene.
- The indications and contraindications of this instructions for use will be considered, in order to minimize error variables.[\*]

11. **Warnings:**

- It is necessary to consider the medical situation of the patient, the physiological and anatomical conditions, and the bone quality and bone quantity in the area where the surgery will be performed. An inadequate technique can incur any of the possible adverse reactions and risks mentioned above.
- The product is perfectly packaged and heat-sealed. Check if the sealed package is damaged or opened because the cleaning conditions of the product could have been altered; in that case, do not use it.
- All products have a batch number to allow their traceability, providing greater security to the intended user, due to with this batch number it is possible to know the entire history of the product from the manufacturing process to the moment of its distribution.
- Any serious incident related to the device must be notified to the manufacturer and to the competent authority of the Member State in which the intended user and/or the patient is established.

12. **Handling:**

- It is important to verify the compatibility of the implant model to be used.
- To place the prosthetic components, compatible dental instruments provided by MPI must be used; to activate the rotation of the power instruments, they must be connected to the handpiece (contra-angle) of the implant motor.

13. **Storage and transport:**
The prosthetic components are wrapped and packaged according to defined and validated protocols that guarantee optimal results. The original packaging has been validated simulating the transport and storage conditions to guarantee the cleanliness of the packaged products, their protection against impacts.
Recommended ideal storage conditions: store in a clean areas, not exposed to direct sunlight.
Incorrect storage during transport can influence the characteristics of the device and lead to its failure.

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