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2. Description and material:

Dental Implants are medical devices (dental), subgingival endosseous implants made of commercially pure Titanium Grade IV cold drawn bars.

The Grade IV Titanium of MPI's Dental Implants meet the chemical composition requirements detailed in ISO 5832-2 standard. Its quantitative information is indicated in the following table:

%C	%N	%O	%Fe	%H
≤ 0.05	≤0.05	≤0.4	≤ 0.5	≤ 0.012

Dental Implants are medical devices class Ib.

3. Intended purpose:

Dental implants intended purpose is to be surgically implanted in the maxilla or jaw, in order to provide fixation to support and retain dental prosthesis to achieve an aesthetic and functional improvement, a better distribution of chewing forces and a decrease in bone loss.

4. Medical purpose / Indication:

Dental implants are indicated for partial or total rehabilitation of edentulous spans using an immediate or delayed loading protocol, according to the intended user's choice, the technique applied and the primary stability achieved.

5. Intended user:

The intended user is a medical professional, responsible for the appropriate surgical procedures and restorative techniques; who must evaluate the suitability of the procedure used, based on their personal medical knowledge, education, training and experience (at least 1 year of experience in dental surgeries with this type of product).

6. Intended patient:

Dental implants are indcated to be placed in: Adult patients (generally older than 17-18 years) with jaw/maxilla growth complete, with unitary spaces without tooth with adjacent and healthy teeth and/or partial edentulism with missing posterior tooth and/or complete edentulism; that cannot tolerate removable prothesis with high aesthetic and/or functional demands, and who are not contraindicated. See "Contraindications" chapter. There is no available data related to special patient populations such as children and during pregnancy or breastfeeding. Therefore, the use of MPI Dental Implants must be according to the criteria of the specialist healthcare professional.

7. Indications for use:

- It is recommended to follow the instructions for use recommended by the manufacturer, as well as to take the necessary protective and clothing measures. Improper use can lead to injure the patient.
- They are intended to be used in a sterile environment.
- The dental implants indicated must be used with the rest of the articles and references of the system to which they belong in order to guarantee a good result.
- Intended users should review the expiration date indicated on the package and check that the package is not damaged or opened to maintain the sterility.
- Dental implants are available in different lengths and diameters in order to be able to rehabilitate different cases according to vertical and horizontal bone availability.
- The implants are delivered packaged with their corresponding cover screw (included in the tube cap). Before placing the implant, the user should verify the appropriate conditions for use them and that the cover screw is suitable for that implant. The colour of the cover screw must match the colour of the line on the box label. Each colour corresponds to a specific diameter, which can be consulted in the catalogue.
- Short® dental implants are recommended for multiple restatations, they are not recommended for immediate loading.
- When it is decided to discard an implant, the intended user must dispose it safely following the protocol for the elimination of medical devices defined in the dental clinic, in accordance with hospital regulations and local existing legislation, to guarantee the safety elimination and avoid infections and contamination.
- Maximum recommended torques depending on the dental implant type (established according to the mechanical tests carried out):

Tradename	Model Number	Maximum recommended torques (N-cm)
MPI	Privilege® EH	110
MPI	Privilege® MT	80
MPI	Excellence®	80
MPI	Short®	80
MPI	All-In®	80

8. Contraindications:

- It is contraindicated to use the dental implants in patients: That suitable for an oral surgical procedure indicated on the package and check that the package is not damaged or opened to maintain the sterility.
- Allergic or hypersensitive to the materials used.
- With periodontal diseases, limited bone availability and density, diabetes and/or smokers and/or bruxism. In case of a patient with bruxism (the involuntary habit of clenching and grinding the teeth without functional purpose), the intended user should assess it in the diagnosis and also in the treatment plan because it should be reflected in theocclusal protection of the implant during bone integration, and in the design of the final prosthesis.
- The placement of dental implants is contraindicated in those patients who cannot undergo surgery for any medical condition.
- The intended user must inform the intended patient of these contraindications.

9. Possible adverse reactions and risks:

- An improper use of the products, an inadequate drilling protocol/surgical technique, can cause bone loss, fenestrations, dehiscences, inflammation, bone and/or nerve damage, discomfort and pain in the patient, and the consequent non-ossseointegration of the implant.
- An inadequate drilling protocol with insufficient irrigation (not cooling the instruments with sterile saline solution) could also cause necrosis in the bone due to thermal effects.
- The inadequate application of higher torque values than the recommended ones, could damage the connection between the implant and the contra-angle; in addition, it could lead to a lack of primary stability, and the bone could

uffer necrosis, which would mean failure in the implant osseointegration.

- The unskilled or untrained user, without planning clinical cases and/or in contraindicated patients, may cause harm to the patient.
- An inappropriate environment where surgery is to be performed/improper storage conditions, can cause contamination of the product, and can lead to infection in the patient and/or non-ossseointegration of the implant.
- To date, no other possible significant adverse effects caused by these products are known.

The intended user must inform the intended patient about the importance of maintaining good dental hygiene, and about these possible adverse reactions and risks.

10. Sterilization:

The indicated dental implants are delivered sterile. They are sterilized by gamma irradiation and they are intended for single use.

11. Precautions:

- Before placing the dental implant, check that it is sterilized (the irradiation label must be red).
- It is recommended that the intended user has received regulated and supervised training by expert professionals on surgical techniques in oral implantology; the intended user should consult the specific surgical and prosthetic protocols; the intended user should receive personalized assistance and the intended user should attend special training before using the products. MPI offers courses for different levels of knowledge and experience.
- The intended user who uses these products must have experience to use them safely and properly according to these instructions for use, in addition the intended user must inform intended patients about contraindications and possible adverse reactions, as well as the importance of maintaining good dental hygiene.
- During the phase of the diagnosis and planning the implant case, the indications and contraindicators of these instructions for use will be considered, the MPI implant system and its prosthetic possibilities should be known, in order to minimize error variables.
- Once the complete case has been planned, and the location, the diameter and the length of the implant have been decided, the type of prosthesis to be used will be designed.
- Subsequently, the patient will undergo surgery, taking into account the design of the incision, the drilling protocol, the placement of the implant and the sterilization conditions of the operative field.

12. Warnings:

- It is necessary to consider the medical situation of the patient, the physiological and anatomical conditions, and the bone quality and bone quantity in the area where the implant is going to be implanted. An inadequate technique can incur any of the possible adverse reactions and risks mentioned above.
- This product is for single use, it cannot be resterilized or re-used. If dental implants are re-used, the sterilization levels offered by MPI are not guaranteed, as they can cause infections and adverse reactions in the patient and therefore the non-ossseointegration of the implant.
- The product is perfectly packaged and heat-sealed. Check if the sealed package is damaged or opened because the sterility and cleaning conditions of the product could have been altered; in that case, do not use it.
- All products have a batch number to allow their traceability, providing greater security to the intended user, due to with this batch number it is possible to know the entire history of the product from the manufacturing process to the moment of its distribution.
- This product is intended to be used in a sterile environment.
- Any serious incident related to the device must be notified to the manufacturer and to the competent authority of the Member State in which the intended user and/or the patient is established.

13. Magnetic resonance safety information:

MPI Dental Implants have not been evaluated for safety and compatibility in the magnetic resonance environment. Scanning a patient with this medical device may not be safe.

14. Handling:

- To place the implants, the implant motor connector must be used; to activate the rotation of this instrument, they must be connected to the handpiece (contra-angle) of the implant motor.
- To accurately measure the drilling mode with the dental drill, it is recommended to verify the drilling depth using calibrated depth gauges.
- To avoid mechanical and thermal risks, use abundant irrigation with sterile saline solution continuously throughout all the implant insertion process.

15. Storage and transport:

The implants are wrapped and packaged according to defined and validated protocols that guarantee optimal results. The original packaging has been validated simulating the transport and storage conditions to guarantee the cleanliness of the packaged products, their protection against impacts, and the sterility of the implants.
Recommended ideal storage in its original package in clean and dry areas avoiding exposed to direct sunlight and for 5 years (expiration period to ensure its sterility). Incorrect storage during transport can influence the characteristics of the device and lead to its failure.

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